



MS[®] Cell/Tissue Culture Plate

Introduction













MS[®] Cell/Tissue Culture Plates are designed under SBS -3D-standard and manufactured of (USP class VI) virgin 100% polystyrene to eliminate all extraneous materials and ensure the integrity of your cells. TC treated or non-treated.



Related types of plates with different well base shape has offered:



Order Information

Type						
Well Format	6	12	24	96	96	96
Base Shape						
Max Volume	17.0 ml	6.2 ml	3.5 ml	380 μ l	380 μ l	380 μ l
Working Volume	2~10 ml	1~3 ml	0.5~1.5 ml	20~340 μ l	20~330 μ l	20~320 μ l
Growth area	9.7 cm ²	3.7 cm ²	2.0 cm ²	0.35 cm ²	—	—
TC Surface Treatment	✓	✓	✓	✓	—	—
γ-irradiated	✓	✓	✓	✓	✓	✓
Cat.-No.	LBCP06S	LBCP12S	LBCP24S	LBCP96S	LBCP96UN	LBCP96VN

Related Products

